

Custom Ptfе Microwave Digestion Vessels For Demanding Applications

Item Number: PL-1003



Introduction

High-purity PTFE microwave digestion vessels for safe, contamination-free sample prep. Ideal for ICP-MS, AAS, and trace analysis. Custom sizes available.

[Learn More](#)

Product Specification	Inner Mouth Diameter (mm)	Body Diameter (mm)	Height with Lid (mm)
PTFE Digestion Vessel 30ml	25	42	105
PTFE Digestion Vessel 50ml	29	47	110
PTFE Digestion Vessel 100ml	50	62	143
PTFE Digestion Vessel 200ml	53	69	163