

Customizable Ptfе Scrapers And Shovels For Demanding Applications

Item Number: PL-1017



Introduction

High-purity PTFE scrapers & shovels for labs, semiconductor & chemical industries. Chemical-resistant, non-stick, durable tools for precise material handling. Custom solutions available.

[Learn More](#)

Type	Width (mm)	Length (mm)	Thickness/Height (mm)	Handle Length (mm)
Scraper	40mm	90mm	5mm	-
Scraper	60mm	120mm	5mm	-
Scraper	100mm	150mm	5mm	-
Scraper	150mm	200mm	5mm	-
Scraper	200mm	200mm	5mm	-
Scraper Board (T-Scraper)	150mm	120mm	5mm	-
Scraper Rod	10mm	120mm	10mm	-
Scraper Strip	15mm	120mm	5mm	-
Shovel (Flat Bottom)	60mm	120mm	30mm	1000mm
Shovel (Flat Bottom)	100mm	120mm	40mm	120mm
Shovel (Flat Bottom)	120mm	150mm	60mm	100mm
Shovel (Flat Bottom)	160mm	200mm	50mm	1000mm
Shovel (Flat Bottom)	170mm	220mm	70mm	115mm
Shovel (Curved Bottom)	110mm	150mm	80mm	100mm
Shovel (Curved Bottom)	110mm	220mm	100mm	120mm
Flat Shovel (Long Handle)	150mm	80mm	10mm	900mm
Flat Shovel (Long Handle)	150mm	100mm	30mm	800mm
Flat Shovel (Long Handle)	200mm	150mm	15mm	500mm