

Dental Porcelain Zirconia Sintering Ceramic Vacuum Press Furnace

Item Number: KT-DP10



Introduction

Precision Vacuum Press Furnace for labs: $\pm 1^{\circ}\text{C}$ accuracy, 1200°C max, customizable solutions. Boost research efficiency today!

[Learn More](#)

Heating speed	200°C/min max. (392°F /minute)
Max temperature	1200°C max (2283°F)
Max temperature holding time	2 hours
Vacuum level	1-101%
Electric power supply	100-120V 50/60HZ 230V 50/60HZ
Rated current	14.5ampere@110V 7.0ampere@230V
Rated output power	max 1500w
Dimension	H27cm x W36cm x D70cm
Firing chamber diameter	9cm(3.5")
Net weight	23Kg
Weight with package	26.5Kg